


BoxPlotter for ChIPseq Regions-of-Interest

Based on the idea in the right picture, we decided to create a more general tool that would allow us to specify sets of regions and compare groups of different samples with each other.



DNase-seq data from 9 male livers (ENCODE) were compared to the STAT5-high/low profiles. **4 livers** showed chromatin accessibility patterns similar to **STAT5-high** livers, while **5 livers** resembled **STAT5-low** livers. These patterns were visualized using read count distributions across the top differential sites.

Algorithm:

1. Calculate FRiP (Fraction of Reads in Peaks)
2. Extract fragment counts from BAM files over target regions for selected samples.
3. Normalize fragment counts using FRiP, generate aggregated and group-wise comparison plots, and compute statistics for group comparisons.

FRiP Normalization:

1. Call MACS2 peaks per sample
2. Merge to union peak set
3. Count fragments per sample in union
4. Normalize to sample-average fragment count

Data transformation during normalization

chr	start	end	width	G241_M01 (Male_2wk_ATAC)			
				RAW COUNTS	RAW_AFTER_WIDTH_NORM	AFTER_FRIP_NORM	FINAL_LOG2_AFTER_RIPPM_NORM
chr1	40293477	40294208	732	706	964.481	395.799	8.629
chr1	40299557	40299948	392	257	655.612	269.047	8.072
chr1	41608251	41608848	598	81	135.452	55.586	5.797
chr1	55136533	55136946	414	384	927.536	380.637	8.572
chr1	82154931	82155388	458	234	510.917	209.667	7.712
chr1	83896719	83897207	489	16	32.720	13.427	3.747
chr1	129796083	129796418	336	158	470.238	192.974	7.592
chr1	133837634	133838008	375	182	485.333	199.169	7.638
chr1	166251542	166251959	418	172	411.483	168.862	7.400
chr1	174860030	174860572	543	173	318.600	130.746	7.031
chr1	176420085	176420403	319	40	125.392	51.458	5.685
chr1	188149410	188149823	414	54	130.435	53.527	5.742
chr10	5728314	5728809	496	147	296.371	121.623	6.926

$1000 * \text{RAW_COUNTS} / \text{WIDTH}$

$\text{RAW_AFTER_WIDTH_NORM} / \text{FRiP}$

- 🏠 Home Directory
- 📁 /project/mzy
- 📁 /projectnb/mzy
- 📁 /restricted/projectnb/mzy
- 📁 /project/wax-arch
- 📁 /projectnb/wax-arch
- 📁 /project/wax-dk
- 📁 /projectnb/wax-dk
- 📁 /project/wax-es
- 📁 /projectnb/wax-es
- 📁 /restricted/project/waxmanlab
- 📁 /restricted/projectnb/waxmanlab

1

2

New Directory

Directory name

your_boxplotter_directory

OK

Cancel

🔍 Open in Terminal ▾

🔄 Refresh

+ New File

📁 New Directory

📶 Upload

⬇ Download

Globus

📄 Copy/Move

🗑 Delete

⬆ / projectnb / wax-es / [Change directory](#)

📄 Copy path


Find your directory and click on it to open

☐ Show Owner/Mode ☐ Show Dotfiles Filter:

Showing 44 rows - 0 rows selected

<input type="checkbox"/>	Type ▲	Name	Size	Modified at
<input type="checkbox"/>	📁	00_shinyapp	-	7/8/2025 11:37:37 AM

Make sure you're in the correct directory and open the Terminal



To clone Boxplotter into the current directory, enter the following command and press Enter:

```
git clone /projectnb/wax-es/WAXMANLAB_SOFT/BoxPlotter_v2 .
```

After cloning the directory you can click to “My interactive sessions” on the top of the page

SCC OnDemand Files ▾ Quotas ▾ Login Nodes ▾ Interactive Apps ▾  My Interactive Sessions 

Interactive Apps

Desktops

Desktop

MATLAB

Mathematica

QGIS

SAS

STATA

Spyder

VirtualGL Desktop

Servers

Jupyter Notebook

MATLAB Server

RStudio Server

Shiny App Server

TensorBoard Server

VS Code Server

Webserver

RStudio Server

This app will launch an RStudio server on a compute node.

Rstudio Version

2024.12.0+467

R Version

4.4.3

Additional modules to load (space separated, optional)

Select Modules

Pre-Launch Command (optional)

Number of hours

12

Number of cores

1

Number of gpus

0

Project

wax-dk

Extra qsub options

☐ I would like to receive an email when the session starts

Launch

Interactive Apps

Desktops

Desktop

MATLAB

Mathematica

QGIS

SAS

STATA

Spyder

VirtualGL Desktop

Servers

Jupyter Notebook

MATLAB Server

RStudio Server

Shiny App Server

RStudio Server (349871)

Host: >_scc-wi4

Created at: 2025-06-11 09:16:54 EDT

Time Remaining: 11 hours and 59 minutes

Session ID: 4eed9c37-ecb2-4287-b0cc-914a422436c3

Rstudio Version: 2024.12.0+467

R Version: R/4.4.3 texlive/2023 pandoc/2.5

Additional modules to load (space separated, optional):

Pre-Launch Command (optional):

Number of hours: 12

Number of cores: 1


Number of gpus: 0


GPU compute capability: 3.5

Project: wax-dk


Extra qsub options:

 Connect to RStudio Server






Click to small arrow to open dialog
and select "Run in Window"



Click "Run App" to start shiny app



BoxPlotter for ChIPseq Regions-of-Interest

Upload Configuration File


Browse...


Choose XLSX file

Select an Excel file (.xlsx or .xls) to use as configuration for the boxplotter job.


Download the example file below to see the required format.

 Download Example File

 Download Instructions



/projectnb/wax-dk/max/src/boxplots_test/rsamtools_atac.Rproj



How to start application

BoxPlotter for ChIPseq Regions-of-Interest

Configuration File

Select XLSX Config file

Browse...

CONFIG_PAIRWISE_


Upload complete

The configuration file must be an XLSX file containing at least three sheets:

- **GROUPS** - defines sample groups
- **COMPARISONS** - specifies comparisons to make
- **REGIONS** - contains regions of interest (can be split into multiple sheets like REGIONS_1, REGIONS_2)

Column order is fixed but names can be arbitrary. Additional columns will be ignored.

Download the example file below to see the required format.

 Download Example File

 Download Results

Once the analysis is complete, download options will appear here.



You can download individual plots or all results as a zip archive.

Once the calculation is complete, a download link will appear.

Analysis Control

 Start Analysis

Example file


 
Starting analysis pipeline... [peak_data]
Calculating FRiP 19 of 21

1


2

3


1. (A vs B) (A_MB_Dynamic_1 (n=834)_top200)




1. (A vs B) (B_MB_Static_Group 1 (n=172))





1. (A vs B) (C_MB Static_Group 2_(n=186))



1. (A vs B) (D_FB_Static_(top n=706)_top200)



A_MB_Dynamic_1 (n=834)_top200





Group name

STAT5_HIGH_Andy

STAT5_LOW_Andy

B_MB_Static_Group 1 (n=172)




Group name

STAT5_HIGH_Andy


STAT5_LOW_Andy

C_MB Static_Group 2_(n=186)

Keep all points for all samples in group for each region



Average of all samples in group for each region




Group name

STAT5_HIGH_Andy


STAT5_LOW_Andy

D_FB Static_(top n=706)_top200

Keep all points for all samples in group for each region



Average of all samples in group for each region




Group name

STAT5_HIGH_Andy

STAT5_LOW_Andy

A_MB_Dynamic_1 (n=834)_top200



B_MB_Static_Group 1 (n=172)

